

Method of reviving from the L-dried ampoule

Prepare three tubes containing a fresh culture medium (5 ml), and label the tubes as No. 1 through No. 3. Suspend the L-dried cells in 0.5 ml of the fresh medium and transfer to tube No. 1. Mix tube No. 1, and transfer 0.5 ml of the suspension to tube No. 2. Then, mix tube No. 2 and transfer 0.5 ml of the suspension to tube No. 3. Incubate all tubes without shaking for 3 to 7 days.

For further information, please contact:  
Takashi ITOH, Ph. D.  
TEL: +81-29-836-9563  
FAX: +81-29-836-9562  
E-mail: [ito@jcm.riken.jp](mailto:ito@jcm.riken.jp)